

The Department of Health has requested modifications to several indicators for nKPI, OSR and PIP-QI reports.

The following table describe the changes that have been implemented in MedicalDirector Insights Version 1.9.0.

Aboriginal Health Program National Key Performance Indicator changes (nKPI)

nKPI Indicator	Summary of changes (as per Specifications for nKPI and OSR, Version 16.0) Specifications for nKPI and OSR — Solving Health
Visit types for nKPI indicators	'Other', 'Telephone' and 'Hospital Consultations' visit types have now been included in Visit Types for nKPI reporting.
PI25 – STI	Addition of local pathology codes and free text search terms to PI25.
PI26- Ear Health Indicator	Inclusion of ear health procedures and conditions as defined in the Ear Health coding framework.
PI11- Smoking status in females	For females, smoking measurements will be considered from the day of the last normal menstrual period (rather the next day which was the previous implementation)

Online Services Reporting (OSR) changes

OSR Indicator	Summary of changes
OSR CS2	Updated CS2 visit type rules to be same as CS1 and CS3

Practice Incentives Program Quality-Improvement (PIP-QI) Indicator changes

Topic/Indicator	Summary of changes
PIPQI JSON file	<ul style="list-style-type: none"> • The 'product_version' field is a dynamic field which reflects the version of MD Insights used to generate the active snapshot • The 'version' field will be 2.0 • The 'author_id' field will automatically display 'Author 1'
QIM3a and b – BMI	<p>We have disregarded BMI measurements and only used height and weight measurements in the calculations for QIM3a and b as per the definition below</p> <p>QIM3a</p> <ul style="list-style-type: none"> • Denominator – Include all regular clients aged 15 years and over. No other exclusions are applied to the denominator. • Numerator – All clients in the denominator if they also had: <ul style="list-style-type: none"> ○ For clients aged 15 years to <25 years old, both height and weight recorded within 12 months. ○ For clients aged 25 years and over, height must be recorded since the client turned 25 and weight must be recorded within 1 year. <p>QIM3b</p> <ul style="list-style-type: none"> • Denominator – 3b denominator should include all clients in the 3a numerator except for clients aged 18 years and over who are shorter than 0.914, or taller than 2.108 meters, or refused measurement. • Numerator – All clients in the 3b denominator who have had their BMI classified under the following categories: <ul style="list-style-type: none"> ○ Obese ○ Overweight ○ Healthy ○ Underweight
QIM8	<p>For the Tobacco Smoking risk factor in the Numerator, we have updated the coding to replace the reference period for 'smoking status recorded within the last 12 months' with 'smoking status ever recorded'</p>